



Instructions for Applying the Dalric Extension #1 and Extension #2

*State of the Art Treatment for
Angular Limb Deformities in Young Foals*



FOOT PREPARATION

- Trim and balance the foot.
- Round the ground surface of the hoof to accommodate the seating of the foot in the cuff.
- Roughen the wall with the smooth side of the rasp. This provides an adequate glue contact surface.
- Clean dirt from all hoof cracks and cover potential abscesses with a small piece of tape or silicone. This prevents sealing bacteria beneath the glue.

CHOOSING THE PROPER SHOE

- Apply the extension lateral (outside) for all toe-in (varus) deformities or medial (inside) for toe-out (valgus) deformities.
- When fitting the shoe, snug the flexible side against the foot. Next, mark the extension and attach the flexible wing to the base of the shoe with the screws provided. Trim excess if necessary.
- Snap the temporary wedge onto the shoe. This keeps the foot snug in the cuff while the glue is curing.

SHOE APPLICATION

- If you are working with an animal that is fractious, you may have your vet administer a sedative.
- Clean the hoof wall with denatured alcohol and ensure the foot is dry.
- Remove the SuperFast™ cartridge and cut the tip off with a pair of nippers. Squeeze a small amount of adhesive from the cartridge to ensure both sides are dispensing equally.
- Apply a mixing tip to the end of the cartridge and twist clockwise to lock it in place.
- Apply SuperFast™ to the felt side of the shoe only.
- Apply the shoe snugly to the clean foot.
- **SuperFast™ sets in 30 seconds, so remember to work quickly.** Total cure time is 2-3 minutes. If you need a longer set time, SuperFast™ can be stored in a refrigerator prior to use.
- To apply a second shoe, simply remove the tip, squeeze a bit of SuperFast™ from the cartridge and apply a new tip.
- When finished, apply a cap to the SuperFast™. The remaining adhesive can be used at another time.

REMOVAL OF THE SHOE

- Use nippers or a rasp to separate the entire cuff from the horn wall and then pull the shoe.

RECOMMENDED SHOETING SCHEDULE

- For the best results, follow this shoeing schedule as the foal's feet continue to grow. Exceeding the recommended shoeing time can alter natural foot development.

Age	Maximum Time for Shoes To Be Worn	Days Following Shoe Removal To Re-evaluate
Birth - 2 Weeks	4 Days	7 Days
2-4 Weeks	7-10 Days	7-10 Days
4-8 Weeks	10-14 Days	7-10 Days
8-12 Weeks	14 Days	7-10 Days
3 Months	21 Days	7-10 Days
8 Months +	28 Days	14 Days

PREVENTION OF SHOE LOSS

- Make sure the shoe fits snugly and the foot is clean prior to application. Keep the foal still during the curing time.
- Do not turn out into mud.
- Do not turn out into large, open fields. Limit exercise to 1-2 acres.
- Use caution when turning out with other horses, as they will step on the extension.

SORENESS FOLLOWING APPLICATION

- Fitting the shoe too tightly will pinch the quarters resulting in slight soreness. If this occurs, pull the shoe and reapply after the foal becomes sound.
- Sealing bacteria in open areas along the white line can cause acute lameness. If this occurs, remove the shoe immediately and notify your veterinarian. The abscess will need to be opened at the white line and cared for appropriately.